

cellSens functions		Dimension	Standard	Entry
Layout	User experience customization	✓	✓	✓
	Overlay multiple images	✓	✓	✓
	Document groups for side-by-side image comparison	✓	✓	✓
	Movie playback	✓	✓	✓
	The view (multiple images in a single data set shown side by side)	✓	✓	✓
View	Slice view for orthogonal plane viewing of 3D or time-lapse data sets	✓	✓	✓
	Voxel view for isosurface and volumetric rendering of 3D and 4D data sets	✓	✓	✓
	Snap/movie acquisition	✓	✓	✓
	Time-lapse at specified interval	✓	✓	✓
	Automated multi-wavelength	✓	✓	✓
Image Acquisition	Z-Stack	✓	✓	✓
	Multi-dimensional (xyzt and wavelength)	✓	✓	✓
	Graphical Experiment Manager	✓	✓	✓
	Manual assisted panoramic imaging (manual MIA)	✓	✓	✓
	Multiposition acquisition and stage navigator	✓	✓	✓
	Automated panoramic imaging (auto MIA, requires motorized stage)	✓	✓	✓
	Instant EFI image (manual or motorized Z)	✓	✓	✓
	Simultaneous multi-color imaging (Image splitter needs)	✓	✓	✓
	Live deblurring	✓	✓	✓
	High Dynamic Range Imaging (HDR)	✓	✓	✓
Image Processing	Multi-well Plate Acquisition	✓	✓	✓
	Geometry/combine/filter processing	✓	✓	✓
	Fluorescence unmixing	✓	✓	✓
	Brightfield unmixing	✓	✓	✓
	Deblurring (No/Nearest Neighbor, Wiener Filter)	✓	✓	✓
Image Analysis	Kymograph	✓	✓	✓
	2D deconvolution	✓	✓	✓
	3D deconvolution (constrained iterative deconvolution)	✓	✓	✓
	Region and line measurements	✓	✓	✓
	Phase analysis	✓	✓	✓
	Object analysis and classification	✓	✓	✓
	Interactive measurement	✓	✓	✓
	Intensity plot over time/z	✓	✓	✓
	Colocalization	✓	✓	✓
	Object Counting (Manual)	✓	✓	✓
Documentation and Collaboration	Online Ratio and Kinetics	✓	✓	✓
	Ratio analysis (off-line)	✓	✓	✓
	FRET analysis	✓	✓	✓
	FRAP analysis	✓	✓	✓
	Automatically compose Word reports	✓	✓	✓
Remoting	Database image and data management solution for microscopy	✓	✓	✓
	Save and load images/documents from Database	✓	✓	✓

Products with confirmed functionality		Dimension	Standard	Entry
Olympus	Camera	DP20 <sup>1</sup> , DP21, DP22, DP25 <sup>2</sup> , DP26, DP27, DP70 <sup>1</sup> , DP71 <sup>3</sup> , DP72 <sup>3</sup> , DP73 <sup>3</sup> , DP80 <sup>3</sup>	✓	✓
	Microscope	BX43, BX53, BX63, BX61, BX61WI, IX83, IX73, IX81, SZX16A	✓	✓
	Peripherals	IX81-ZDC, IX81-ZDC2, IX3-ZDC, IX3-ZDC2	✓	✓
	Motorized XY stage	BX-DSU, IX3-DSU, IX2-DSU, U-CBF	✓	✓
Olympus Soft Imaging Solutions	Camera	BX3-SSU, IX3-SSU	✓	✓
	Peripherals	CC12, F-View II, Colorview I, Colorview II, Colorview III, Colorview IIIu, XM10, XC10, XC30, XC50, UC30, UC50, SC20, SC30, SC100	✓	✓
Hamamatsu	Camera	cellTRF (multi-line, single line), MT20, USB-ODa converter, Real Time Controller (U-RTC and U-RTCE), U-FCB, U-STC	✓	✓
	Image Splitter	ORCA R2, ORCA 03, ORCA 05, ORCA ERG, ORCA-Flash 2.8, ImagEM, ImagEMX2, ORCA-Flash 4.0 V2, ORCA-Flash 4.0 LT	✓	✓
Q-Imaging	Camera	W-View Gemini	✓	✓
	Image Splitter	MicroPublisher 3.3 RTV, MicroPublisher 5 RTV	✓	✓
Photometrics	Camera	Monochrome: EXi Blue/Aqua, QiClick, Retiga EXi/SPV/2000R/2000RV/4000R/4000RV/6000	✓	✓
	Image Splitter	Color: iExi Aqua, OptiMOS, Rolera Thunder, CoolSNAP HQ2, Evolve 512 Delta	✓	✓
Andor	Camera	Dual View DV2 /QuadView QV2	✓	✓
	Image Splitter	MicroPublisher 3.3 RTV, MicroPublisher 5 RTV	✓	✓
Jenoptik	Camera	MicroPublisher 3.3 RTV, MicroPublisher 5 RTV	✓	✓
	Image Splitter	MicroPublisher 3.3 RTV, MicroPublisher 5 RTV	✓	✓
Vincent Associates	Camera	MicroPublisher 3.3 RTV, MicroPublisher 5 RTV	✓	✓
	Image Splitter	MicroPublisher 3.3 RTV, MicroPublisher 5 RTV	✓	✓
CoolLED	Light Source	MicroPublisher 3.3 RTV, MicroPublisher 5 RTV	✓	✓
	Image Splitter	MicroPublisher 3.3 RTV, MicroPublisher 5 RTV	✓	✓
Excelitas	Light Source	MicroPublisher 3.3 RTV, MicroPublisher 5 RTV	✓	✓
	Image Splitter	MicroPublisher 3.3 RTV, MicroPublisher 5 RTV	✓	✓
Sutter	Light Source	MicroPublisher 3.3 RTV, MicroPublisher 5 RTV	✓	✓
	Image Splitter	MicroPublisher 3.3 RTV, MicroPublisher 5 RTV	✓	✓
Prior	Shutter, FW	MicroPublisher 3.3 RTV, MicroPublisher 5 RTV	✓	✓
	Motorized XY stage	MicroPublisher 3.3 RTV, MicroPublisher 5 RTV	✓	✓
Ludl	Shutter, FW, Z-drive	MicroPublisher 3.3 RTV, MicroPublisher 5 RTV	✓	✓
	Piezo Z (control via Real Time Controller)	MicroPublisher 3.3 RTV, MicroPublisher 5 RTV	✓	✓
Objective Imaging	Motorized XY stage	MicroPublisher 3.3 RTV, MicroPublisher 5 RTV	✓	✓
	Shutter, FW, Z-drive	MicroPublisher 3.3 RTV, MicroPublisher 5 RTV	✓	✓
Märzhäuser	Motorized XY stage controller	MicroPublisher 3.3 RTV, MicroPublisher 5 RTV	✓	✓
	Z-drive controller	MicroPublisher 3.3 RTV, MicroPublisher 5 RTV	✓	✓
Physik Instrumente	Motorized XY stage	MicroPublisher 3.3 RTV, MicroPublisher 5 RTV	✓	✓
	Piezo Z (control via Real Time Controller)	MicroPublisher 3.3 RTV, MicroPublisher 5 RTV	✓	✓
Applied Scientific Instrumentation	Motorized XY stage	MicroPublisher 3.3 RTV, MicroPublisher 5 RTV	✓	✓
	Piezo Z (control via Real Time Controller)	MicroPublisher 3.3 RTV, MicroPublisher 5 RTV	✓	✓
National Instruments	Digital TTL device	MicroPublisher 3.3 RTV, MicroPublisher 5 RTV	✓	✓
	CSU	MicroPublisher 3.3 RTV, MicroPublisher 5 RTV	✓	✓
Yokogawa	Digital TTL device	MicroPublisher 3.3 RTV, MicroPublisher 5 RTV	✓	✓
	CSU	MicroPublisher 3.3 RTV, MicroPublisher 5 RTV	✓	✓

Compatible image formats	
Read and write	JPEG, JPEG2000, TIFF, BMP, AVI, PNG, VSI (Virtual slide image),
Read only	GIF, PSD (Adobe Photoshop), TIFF (DP-BSW, FSX100, MetaMorph), OIF/OIB/OIR (Fluoview format), Cell, STK (MetaMorph), MRC (Medical Research Council)

Recommended system requirements	
OS	Microsoft Windows 8.1 Pro (32-bit/64-bit), Microsoft Windows 8 (32-bit/64-bit) Pro, Microsoft Windows 7 (32-bit/64-bit) Ultimate with SP1, Microsoft Windows 7 (32-bit/64-bit) Professional with SP1
OS Language	English, Simplified Chinese, Japanese, German, Russian (only for Entry and Standard) and all others with English like alphabet
CPU	Intel Core i5, Intel Core i7, Intel Xeon Recommended for high speed image acquisition: QuadCore
RAM	4 GB Recommended for high speed image acquisition: 8GB or more only on Windows7 64-bit operating system
Graphic card	1280x1024 (min. 1024 x768) monitor resolution with 32-bit-video card with separate graphics memory (no integrated graphics processor with shared memory)
Port	USB 2.0 port to connect devices to the system FireWire A to connect devices to the system (BX61, IX81, SZX2-MDCU, IX3-DSU etc...) Serial (RS232) to connect devices to the system (BX61, IX81, SZX2-MDCU etc...) Additional PC/PCIe slots as necessary to connect third party peripherals (principally third party cameras) with propriety interface cards
HDD	1 GB for installation Performance of hard disk is a limiting factor for image acquisition speed Recommended for high speed image acquisition: Solid State Drive (SSD)
Drive	DVD drive (Read: DVD-R DL)
Web Browser	Recommended for Windows 7: Microsoft Internet Explorer 8.0, 9.0, Recommended for Windows 8: Microsoft Internet Explorer 10, Recommended for Windows 8.1: Microsoft Internet Explorer 11

Image data courtesy of:  
Hiroo Ueno, Ph.D.  
Department of Stem Cell Pathology, Kansai Medical University  
(cover page)

[www.olympus-lifescience.com](http://www.olympus-lifescience.com)



For enquiries - contact  
[www.olympus-lifescience.com/contact-us](http://www.olympus-lifescience.com/contact-us)

**OLYMPUS CORPORATION**  
Shinjuku Monolith, 2-3-1 Nishi-Shinjuku, Shinjuku-ku, Tokyo 163-0914, Japan  
**OLYMPUS EUROPA SE & CO. KG**  
Wendenstrasse 14-18, 20097 Hamburg, Germany  
**OLYMPUS SCIENTIFIC SOLUTIONS AMERICAS CORP.**  
48 Woerd Avenue, Waltham, MA 02453, U.S.A.  
**OLYMPUS SINGAPORE PTE LTD.**  
491B River Valley Road, #12-01/04 Valley Point Office Tower, Singapore 248373  
**OLYMPUS MEDICAL SYSTEMS INDIA PRIVATE LIMITED.**  
102-B, First Floor, Time Tower, M.G. Road, Gurgaon 122001, Haryana, INDIA


**OLYMPUS LATIN AMERICA, INC.**  
5301 Blue Lagoon Drive, Suite 290 Miami, FL 33126, U.S.A.  
**OLYMPUS (CHINA) CO., LTD.**  
A8F, Ping An International Financial Center, No. 1-3, Xinyuan South Road, Chaoyang District, Beijing, 100027 P.R.C.  
**OLYMPUS KOREA CO., LTD.**  
8F Olympus Tower, 446 Bongseunsa-ro, Gangnam-gu, Seoul, 135-509 Korea  
**OLYMPUS AUSTRALIA PTY. LTD.**  
3 Acacia Place, Notting Hill VIC 3168, Australia

Printed in Japan M1718E-092015



Your Vision, Our Future

Seamless Workflow. Intuitive



Imaging Software

cellSens

Not for clinical diagnostic use.



# Add Simplicity to Experiment Design... Leave More Time for Research


Olympus cellSens gives you a simpler way to work.

Enjoy full control over the user interface, with functions that are where you want them, when you need them.

Seamless operation, from image capture to report creation means more results with less effort.

Spend less time with your software. Have more time for research.


## Imaging



### Image Capture

Capture multi-color, time lapse, and z-stack images with ease. Just select the appropriate capture button, add relevant parameters, and click "Start". The Process Manager or Experiment Manager make it easy to capture multidimensional image.


## Processing



### Viewing and Processing

Automatically view your data in the colors and layout you choose. Take advantage of an array of advanced image processing functions, such as stitching, extended focus, deconvolution, and unmixing.


## Analyzing



### Measurement and Analysis


Make measurements using an intuitive interface. cellSens offers region of interest, phase analysis, and cell count capability. Export raw measurement data to MS Excel or a cellSens workbook with a single click.

## Reporting



### Collaboration and Communication

Actively collaborate with colleagues and coworkers with special tools including Database and Reporting functions. These functions make it simple to manage, share, and distribute your own image and data reports.



## Microscopy Research With a Personal Touch

Olympus microscopes enable new techniques and push the boundaries of resolution at all magnifications. Olympus cellSens software improves productivity with efficient acquisition work-flow, image processing capability and analytical strength. Centered around the needs demanding customers, cellSens is flexible, easily personalized, and designed to fit as application requirements evolve.



# Reduce Clutter and Confusion by Displaying Only the Tools and Windows You Need

### It's Time to Get Personal

Olympus has been at the forefront of microscopy for over 90 years and has developed microscopes and systems for a broad spectrum of applications. As a result, we know that each researcher has individual requirements that can't all be met by fixed solutions. The cellSens software family consists of three packages, all featuring a peerless user-definable interface. As a result, each user can define what they want cellSens to show them within the defined work areas.

### Dynamic Interface


Creating an efficient workflow requires careful definition of the tasks and tools at each stage. With the cellSens platform's dynamic GUI, the same is true—the tools you need for each stage are clearly available, without clutter or the need to search. Olympus has created a number of interface layouts, which are developed with capabilities appropriate to the users needs.

- **Acquisition Layout**—for selecting between different acquisition processes and adjusting the camera settings.
- **Processing Layout**—for post-acquisition functions such as image processing, execution of measurements, collection of data, presentation of resulting statistics.
- **Count & Measure Layout**—for manual and automated measurement and object counting.
- **Reporting Layout**—for generating reports to document and share results.
- **Create Layout**—a user can define his or her own layout in various arrangements.




### Camera Control Panel

The most important microscope component that requires software control when imaging is the digital camera. Modern cameras feature a number of functions that can be changed to enhance or perfect an image; for example, exposure time and pixel binning. The cellSens Entry and Standard packages control such features on all Olympus digital microscopes and cameras. The Dimension package, in addition, controls such features on high-end research cameras as well. As a result, scientists can maximize the quality of their images.



### Dark Application Skin

The Dark Application Skin reduces computer monitor-generated ambient light and allows cellSens users to adapt to darkened environments; icon contrast remains high for easy recognition and quick selection.




Dark skin




### Arrange Windows as You Like


Organize the tools and windows for the job at hand to create a functional layout that works best for you.



Full screen



Floating panels



Docked panels

Need Help? Online Help is Just a Click Away

Display Only Those Functions You Need on the Toolbar

Common Functions can be Grouped in a Single Tab

All necessary functions are placed where you want, when you need them. Layout tabs allow easy selection of functions according to your workflow. For instance, display camera control features in your Acquisition layout, and then remove them from view when you switch to the Processing layout.

Display or Hide Windows as You Require, or Use Auto-hide for Clean Operation

Graphical Experiment Manager (GEM)

GEM enables the design of complex experiments by simply dragging and dropping icons onto the canvas.

Create Flexible Workflow Toolbars for Repetitive Operations

cellSens lets you create custom toolbars for your most frequently used functions and then save them to the My Functions window. Custom buttons are also easy to use, with convenient tab access that further enhances workflow efficiency. Furthermore, appearance of each toolbar in the My Functions window can be customized by choosing an icon and/or text from Button Appearance window.

Functional Panels are Grouped in Tabs for Easier Access


# Solutions to Empower Your Individual Research

## What Scientific Researchers Wanted

## Our Solutions


Complex experimental procedure with flexible design

**Graphical Experiment Manager (GEM)**  
Experiments can be freely designed simply by connecting the various commands. Furthermore, image acquisition is available for up to 6 dimensions (XYZTλ multipoint).




Flexible well plate image capture

**Well Plate Navigator**  
Capture well plate samples automatically by using the well plate navigator in combination with the motorized stage. There is also enhanced flexibility to allow multiple experiments to be executed within a single well plate.




Intensity analysis

**Intensity Analysis**  
Visualize changes in intensity over time, and save this information for later analysis. Ratio Analysis function allows calibration, display and analysis of live/stored data reflecting changes in the intensity ratio between two acquisition channels.



Improved image detail


**Deconvolution**  
Choose between 2D (included) and 3D (optional) blind deconvolution. This proprietary and highly efficient post-processing tool for both CCD and Confocal imaging enhances the ability to differentiate between imaged objects.



Kei Ito, Ph. D.  
Institute of Molecular and Cellular Biosciences, University of Tokyo

Unified task order management

**Macro Manager**  
Perform tasks, from imaging to processing and analysis, as a single composite procedure. Batch processing is also available, enabling multiple images to be subjected to preferred processes as a continuous series for a significant improvement in workflow efficiency.




## What Medical Researchers Wanted

## Our Solutions

Retention of intact observed images

**Live/Snapshot Function with White Balance Adjustment**  
Simply align the focus and select the appropriate white balance to capture images with true-to-life quality.




Observe large sample at once

**Panoramic Imaging**  
Create clear and seamless wide area images by automatic correction of mismatching between each images, even when using the manual stage. A fully functional wide-area focus map enables improved clarity in panoramic imaging.




Simultaneously monitoring of multiple images

**Image Comparison (simultaneous image windows)**  
Display images side by side for accurate comparison, with simultaneous zooming and movement.




Cell counting by hands

**Object Counting**  
Perform manual counts with self-set classes. Counts and proportions can then be undertaken for each class through simple mouse operation.



Nuclei counting with variable thresholding

**Particle Analysis**  
Set threshold levels for nuclei counts, or calculate parameters such as tissue slice total area and area ratios.





# An Array of Easy-to-Use Functions


## Turn Research Finding into Compelling Presentations

### Image Capture


#### Graphical Experiment Manager (GEM)

Dimension

A break from the usual complex panel-based interfaces, the Graphical Experiment Manager (GEM) uses a flexible drag-and-drop interface to build simple or complex experiments within the cellSens workspace. Actions can be combined within specialized frames to dictate the desired order and priority of automation and imaging. Easily acquire multichannel imaging, Z-stacks, or time-lapse acquisitions across one or more sample positions. GEM permits user interaction with the system during automation to address unforeseen changes in the sample, save time, and prevent repeat effort.



Capture command




Loop command

#### Well Plate Navigator

Dimension + Multiposition + Well Plate Navigator


The Well Plate Navigator Solution automatically scans and acquires images from standard and customized well plate formats. All acquired images, sample positions, and user comments can be saved into a structured database for rapid centralized access. Unique imaging settings can be flexibly applied to individual wells, columns, and rows. Additionally, the Well Plate Navigator supports the execution of multiple experiments within a single well plate to support more complex experiments.



#### Capture Multidimensional Images

Dimension + Multiposition

In combination with a motorized microscope, the Process Manager makes it easy to capture multi-colored and multidimensional images. With the optional Multi-position Solution, it is also possible to capture multi-point and large area images automatically.



Process Manager Setting


Multi-color images

#### Panoramic Imaging

Dimension + Multiposition

Standard + Manual Process

The manual multiple image alignment function creates a single panoramic image as the specimen is scanned. Wide area imaging using a motorized stage can also be fully automated with cellSens Dimension and the optional Multi-position Solution. In combination with a motorized z-focus, this function captures images that are auto-corrected for sample distortion and tilting.



#### Extended Focus Imaging

Dimension

Standard + Manual Process


cellSens will create a single in-focus image from successive image planes as the focus knob is turned using the Extended Focus Imaging (EFI) function, or motorized focus drive to fully automated EFI acquisition. EFI composites images can also be created directly from previously captured Z-stacks.

### Viewing and Processing

#### Unmixing

Dimension

The linear unmixing algorithm in cellSens Dimension results in crosstalk-free fluorescent images to address the challenge of fluorochromes with overlapping emission spectra - such as GFP and YFP. cellSens linear unmixing can also help to separate background autofluorescence from fluorescence signal.




Before unmixing

After unmixing

#### Best Focus Extraction

Dimension

Extract the best focus from images, including z-stack, time-lapse images. This function is effective in creating T-series images with the best focus possible, even when working with defocused time-lapse images.



#### High Dynamic Range Imaging (HDRI)


Dimension

By automatically capturing many images at different exposures the HDRI function creates a final image with a much greater dynamic range, where low intensity signals are clearly visible without overexposing the bright areas of the sample.

#### Deconvolution

Dimension + CI Deconvolution

The optional CI Deconvolution Solution employs the latest in Constrained Iterative Deconvolution algorithms to produce improved resolution, contrast and dynamic range, with industry-leading speed. Each cellSens Dimension license includes Live 2D deblurring for preview and acquisition to enable better contrast within thick specimens. cellSens comes complete with the most widely requested deblurring techniques such as 2D deconvolution, Nearest neighbor, Wiener filter.



Original image

Deconvolved image

Cell line: Human cervical cancer cell line HeLa cell  
Immunostaining: Hec1 staining (green, Alexa Fluor-488),  $\alpha$ -tubulin staining (red, Alexa Fluor-568), DAPI staining (blue)  
Mitotic HeLa cell derived from human cervical cancer.  
Mitotic spindle and kinetochores are stained with anti- $\alpha$ -tubulin (red) and anti-Hec1 (green) antibodies, respectively. Chromosomes interact with microtubules constituting mitotic spindle via kinetochores, protein structure assembled on centromere region of chromosomes.

**Image data courtesy of:**  
Department of molecular oncology, Institute of Development, Aging and Cancer, Tohoku university  
Masanori Ikeda and Kozo Tanaka

#### Solution

Each cellSens Package can be expanded towards a specific application by using optional "Solutions"

Dimension	available solutions:	
CI Deconvolution	Multiposition	Well Plate Navigator
Count & Measure	Ratio/FRET	Database Core
Database Client	NetCam	Photo Manipulation
Standard	available solutions:	
Multichannel Acquisition	Manual Process	Database Core
Database Client	NetCam	
Entry	available solution:	
Database Client		

# An Array of Easy-to-Use Functions


## Turn Research Finding into Compelling Presentations

### Measurement and Analysis

#### Manual Measurement

- Dimension
- Standard
- Entry



Distances between points, areas, intensity measurements and morphological parameters become accessible using the cellSens measurement tools. Measurement data are saved as an image layer that can be exported to Microsoft Excel and cellSens workbook formats, or viewed using OlyVIA, the free image viewer software package.



#### Intensity Analysis

- Dimension


Graphically depict intensity and channel ratios, and export values to Excel or WorkBook by simply setting the region of interest (ROI). The ROI can be moved to capture measurements in line with cell movements. Convert variations of intensity to hue and brightness using Intensity Modulated Display (IMD) to visually enhance fine image structure often found within ratio or FRET images. The Ratio/FRET Solution is used for analysis and display of real-time ratiometric imaging and data. FRET analysis of both sensitized emission and acceptor photo-bleaching is also supported within this user friendly work-flow. The Photo-Manipulation Solution can be used for the curve-fitting analysis of FRAP images.



#### Automatic Object Measurement and Classification

- Dimension + Count & Measure


The Count & Measure Solution adds efficient and precise object detection for automated nuclei counting and classification. This solution expands the set of manual measurements in cellSens. Perform automatic object measurement and classification easily, using an interactive interface where recognized objects are always linked with their measurements.



#### Manual Count

- Dimension
- Standard

Perform manual counts with self-set classes. Counts and proportions can then be undertaken for each class through simple mouse operation.




### Collaboration and Communication

#### Reporting

- Dimension

Easily drag-and-drop image property data, measurement data and user-customized fields into a report template using the convenience of a built-in reporting tool to produce Microsoft Word\* format reports. Collaborate with colleagues, and communicate results, quickly and easily.


*\*Requires Microsoft Word version 2007 or later*



#### Remote Live Image

- Dimension + NetCam
- Standard + NetCam


The cellSens NetCam Solution facilitates the transfer of live or static imaging over a network for teaching, mentoring or supervision.



#### Database

- Dimension + Database Core or Database Client
- Standard + Database Core or Database Client
- Entry + Database Client

Database Core Solution enables the creation of network shared, user-definable databases, with full access controls. The database stores images and all associated image properties, user comments and any other related files that a user wishes to include. The interactive query tool makes it easy to find the data, and provides automatic previews of each queried image. Read and write to a shared database conveniently from multiple different stations with the Database Client Solution.



#### Database+well Navigator

- Dimension + Multiposition + Well Plate Navigator
- Database Core or Database Client

In combination with the Well Navigator Solution, the Database Solution greatly improves the efficiency of viewing and analysis of well plate images with a large amount of data. By clicking on icons for image information such as the date, file name, or well plate number, any selection of captured images can be viewed for further analysis. This solution also allows viewing of captured images and continuous analysis of selected images (the Batch Macro function) via the well plate GUI.

